


Pressure Cycling Technology (PCT) Reduces Effects of Inhibitors of the PCR

Pamela L. Marshall¹ · Jonathan L. King¹ · Nathan P. Lawrence² · Alexander Lazarev² · Vera S. Gross² · Bruce Budowle¹

¹Institute of Applied Genetics, Department of Forensic and Investigative Genetics, UNT Health Science Center, 3500 Camp Bowie Blvd., Fort Worth, TX 76107, USA
²Pressure BioSciences Incorporated, 14 Norfolk Avenue, South Easton, MA 02375, USA



Abstract

A common problem in the analysis of forensic human DNA evidence, or for that matter any nucleic acid analysis, is the presence of contaminants or inhibitors. Contaminants may co-purify with the DNA and inhibit downstream PCR, or they may present samples effectively as containing fewer templates than exist in the PCR, even when the actual amount of DNA is adequate. Typically, these challenged samples exhibit allele imbalance, allele dropout and sequence specific inhibition which can lead to interpretational difficulties. Lessening the effects of inhibitors may increase the effective yield of challenged low template copy samples. High pressure may alter some inhibitors and render them less effective at reducing the yield of PCR products. In an attempt to enhance the amplicon yield of inhibited DNA samples, pressure cycling technology (PCT) was applied to DNA exposed to various concentrations of hematin (0, 1.25, 2.5, 5, and 7 μ M) and humic acid (0, 1.25, 2.5, 5, and 7 ng/ μ l). The effect of high pressure on the inhibitors and, subsequently, the PCR process was assessed by measuring DNA quantity by qPCR and evaluating STR typing results. The results support that pressure cycling technology reduces inhibitory effects and, thus, in effect enhances yield of amplified products of both hematin and humic acid contaminated samples. Based on the results obtained in this study, this method can improve the ability to type challenged or inhibited DNA samples.

Materials and Methods

DNA
 Experiments were performed using the Quantifiler[®] Human DNA Quantification Kit (Life Technologies) DNA Standard (Raji cell line; 200 ng/ μ l purified DNA). The standard was diluted to a final concentration of 1 ng/ μ l for each experimental sample.

PRESSURE CYCLING TECHNOLOGY
 Samples undergoing PCT were placed in either FT500-ND PULSE Tubes or MicroTubes. Samples were transferred to the Barocycler[®] NEP3229 and subjected to 30 cycles of alternating pressures consisting of 35kpsi for 20 seconds and ambient pressure for 10 seconds. Non-pressure treated controls (NPC) controls also were prepared.

QUANTIFICATION
 Quantity of DNA was determined using the Quantifiler[®] Human DNA Quantification Kit on the ABI 7500 Real-Time PCR System (Life Technologies).

BONE SAMPLES
 Five human bones were obtained and prepared for DNA extraction. The bones were cleaned and air dried overnight. Then the bones were crushed to powder using a 6750 Freezer/Mill (SPEX SamplePrep L.L.C.), filled with liquid nitrogen, using a protocol of a 10 minute re-chill followed by 5 minutes of grind time at 15 impacts per second. Approximately 0.2 g of bone powder was placed in a PULSE[™] Tube. 1 ml of extraction buffer containing 0.5 M EDTA pH 8.0 (Invitrogen Corporation), 1% sodium lauryl sarcosinate (sarkosyl, n-lauryl sarcosine; Sigma), and 100 μ g/ml Proteinase K (Invitrogen Corporation) were added to each sample and vortexed. Both pressure treated and NPC samples were placed into PULSE Tubes. For three of the bones, the NPC sample contained 0.5 g bone powder. The bones were incubated at 56°C with constant agitation for either two hours or overnight. Following incubation, samples were subjected to either PCT or no pressure. All bone samples were processed using a hi-flow silica column extraction.

AMPLIFICATION AND STR TYPING
 Approximately one nanogram of DNA for each extract was amplified using the reagents contained in the AmpFSTR[®] Identifier[®] PCR Amplification Kit (Life Technologies) according to the manufacturer's recommendations. A subset of experiments was performed using the reagents contained in AmpFSTR[®] Identifier[®] Plus PCR Amplification Kit (Life Technologies) according to the manufacturer's recommendations. PCR products were separated and detected on an ABI 3130xl Genetic Analyzer (Life Technologies) following the manufacturer's recommendations. The detection and interpretation thresholds both were set at 50 relative fluorescence units (RFU).

Results




Figure 1 Effect of PCT on Hematin Inhibition
 Various concentrations of hematin (0, 1.5, 2.5, 5, and 7 μ M) were added to 100 μ l DNA (1 ng/ μ l), placed in MicroTubes and were subjected to either barocycling (30 Cycles [20s at 35kpsi and 10s at ambient psi]) or no pressure. 1 μ l of each sample subsequently was amplified and typed for STRs. Samples were run in replicates of five. Non-pressure cycled samples are on the left with pressure cycled samples on the right and concentration of hematin is as follows: (a) 0 μ M, (b) 1.5 μ M, (c) 2.5 μ M, (d) 5 μ M, and (e) 7 μ M. Pressure samples were compared with non-pressured controls.




Figure 3 Effect of Ethyl Alcohol Addition on Hematin and Humic Acid Inhibition
Top Panel: 2.5 μ M/ μ l of hematin were added to 100 μ l DNA (1 ng/ μ l), placed in MicroTubes and were subjected either to barocycling (30 Cycles [20s at 35kpsi and 10s at ambient psi]) or no pressure. Comparison samples for pressure and no pressure groups were created with the addition of EtOH at a final concentration of 10%. 1 μ l of each sample subsequently was amplified and typed for STRs. Samples were run in replicates of three. Non-pressure cycled samples containing 1 ng/ μ l control DNA and 2.5 μ M/ μ l hematin, without (a) and with (b) addition of 10% Ethyl Alcohol. PCT samples containing 1 ng/ μ l DNA and 2.5 μ M/ μ l hematin, without (c) and with (d) addition of 10% Ethyl Alcohol are displayed.
Bottom Panel: 2.5 ng/ μ l of humic acid were added to 100 μ l DNA (1 ng/ μ l), placed in MicroTubes and were subjected to either barocycling (30 Cycles [20s at 35kpsi and 10s at ambient psi]) or no pressure. Comparison samples for pressure and no pressure groups were created with the addition of EtOH at a final concentration of 10%. 1 μ l of each sample subsequently was amplified and typed for STRs. Samples were run in replicates of three. Non-pressure cycled samples containing 1 ng/ μ l DNA and 2.5 ng/ μ l humic acid, without (a) and with (b) addition of 10% Ethyl Alcohol are displayed. Pressure cycled samples containing 1 ng/ μ l DNA and 2.5 ng/ μ l humic acid exposed to PCT without (c) and with (d) addition of 10% Ethyl Alcohol are displayed.




Figure 2 Effect of PCT on Humic Acid Inhibition
 Various concentrations of humic acid (0, 1.5, 2.5, and 5 ng/ μ l) were added to 100 μ l DNA (1 ng/ μ l), placed in MicroTubes and were subjected to either barocycling (30 Cycles [20s at 35kpsi and 10s at ambient psi]) or no pressure. 1 μ l of each sample subsequently was amplified and typed for STRs. Samples were run in replicates of five. Non-pressure cycled samples are on the left with pressure cycled samples on the right and concentrations of humic acid are as follows: (a) 0 ng/ μ l, (b) 1.5 ng/ μ l, (c) 2.5 ng/ μ l, and (d) 5 ng/ μ l. Pressure samples were compared with non-pressured controls.




Figure 4 Effect of PCT on Human Bone Processing
 Electropherograms representing a human bone sample processed with and without PCT. Two different quantities of bone powder were used for the NPC sample (0.2 g and 0.5 g) while 0.2 g of bone powder was used for the pressure treated sample. The bones were incubated in PULSE Tubes with extraction buffer at 56°C overnight. Following incubation, the pressure treated sample was subjected to barocycling (30 Cycles [20s at 35kpsi and 10s at ambient psi]). Subsequently, DNA was extracted from all samples using hi-flow silica column extraction. The first elution of each sample was amplified (all samples yielded less than 0.1 ng/ μ l, 10 μ l were amplified) and typed for STRs using the AmpFSTR[®] Identifier[®] Plus PCR Amplification Kit.

Conclusions and Future Directions

This study shows the potential of enhanced PCR efficiency for samples containing an inhibitor when PCT treated compared with those samples not exposed to PCT. This research study suggests that PCT potentially has applications for forensic DNA analysis of certain challenged forensic DNA samples by reducing the effects of inhibitors known to be present in some bone samples. Future research will focus on elucidating the mechanism(s) that overcomes the effect of inhibition.

References

Alaeddini R. Forensic implications of PCR inhibition-A review. *Forensic Sci Int Genet.* 2011 Sep 12.
 Eiler KD, Foran DR. Polymerase resistance to polymerase chain reaction inhibitors in bone*. *J Forensic Sci.* 2009;54(5):1001-7.
 Ramakers C, Ruijter JM, Deprez RH, Moorman AF. Assumption-free analysis of quantitative real-time polymerase chain reaction (PCR) data. *Neurosci Lett.* 2003 Mar 13;339(1):62-6.
 Funes-Huacca ME, Opel K, Thompson R, McCord BR. A comparison of the effects of PCR inhibition in quantitative PCR and forensic STR analysis. *Electrophoresis.* 2011;32(9):1084-9.
 Kontanis EJ, Reed FA. Evaluation of real-time PCR amplification efficiencies to detect PCR inhibitors. *J Forensic Sci.* 2006;51(4):795-804.
 Huggert JF, Novak T, Garson JA, Green C, Morris-Jones SD, Miller RF, et al. Differential susceptibility of PCR reactions to inhibitors: An important and unrecognized phenomenon. *BMC Res Notes.* 2008 Aug 28;1:70.
 Budowle B, Eisenberg AJ, van Daal A. Validity of low copy number typing and applications to forensic science. *Croat Med J.* 2009 Jun;50(3):207-17.
 Irwin JA, Leney MD, Loreille O, Barritt SM, Christensen AF, Holland TD, et al. Application of low copy number STR typing to the identification of aged, degraded skeletal remains. *J Forensic Sci.* 2007 Nov;52(6):1322-7.
 Gross V, Carlson G, Kwan AT, Smejkal G, Freeman E, Ivanov AR, et al. Tissue fractionation by hydrostatic pressure cycling technology: The unified sample preparation technique for systems biology studies. *J Biomol Tech.* 2008 Jul;19(3):189-99.
 Smejkal GB, Robinson MH, Lawrence NP, Tao F, Saravis CA, Schumacher RT. Increased protein yields from *Escherichia coli* using pressure-cycling technology. *J Biomol Tech.* 2006 Apr;17(2):173-5.
 Smejkal GB, Witzmann FA, Ringham H, Small D, Chase SF, Behnke J, et al. Sample preparation for two-dimensional gel electrophoresis using pressure cycling technology. *Anal Biochem.* 2007 Apr 15;363(2):309-11.

Acknowledgements

This project was supported by Award No. 2009-DN-BX-K188, awarded by the National Institute of Justice, Office of Justice Programs, U.S. Department of Justice. The opinions, findings, and conclusions or recommendations expressed in this publication/program/exhibition are those of the author(s) and do not necessarily reflect those of the Department of Justice.

Poster Presentation at International Symposium for Human Identification, October 2012, Nashville, TN

Contact Information

Contact Pam Marshall for information regarding the content or a reprint of this poster at: pamela_marshall@untthsc.edu

Article Information

This research was accepted for publication in the International Journal of Legal Medicine. The complete article can be accessed via the following DOI: <http://dx.doi.org/10.1007/s00414-012-0770-y>