


# Lysis Rate of *Escherichia coli* With Varying Pressure and Cell Count

Lysis of *Escherichia coli* at Different Cell Counts and at Increasing Pressures




**Culture:** A single colony of *Escherichia coli* ATCC 8739 was used to inoculate 10 mL sterile LB (Lysogeny Broth: Tryptone 10 g/L, Sodium Chloride 10 g/L, Yeast Extract 5 g/L) and was grown overnight at 37°C with shaking at approximately 180 rpm. This starter culture was used in a ratio of 1% to inoculate 200 mL sterile LB in 1 L flasks. These were then grown under the same conditions until an Optical Density of approximately 0.6 was reached ( $10^8$  cells/mL). At this point the culture was split; half was refrigerated and half was grown further until Optical Density of 1.2 was reached.

**Disruption:** The resulting culture was passed in 30 mL aliquots through a Constant Systems Ltd. 'Z+' 1.1 kW Continuous Flow Cell Disruptor at pressures indicated. The machine was rinsed with 30 mL de-ionised water between each use.

**Counting:** After being passed through the machines, 10 µL of lysate was mixed 1:1 with the viability stain Trypan Blue. Live cells were counted using a hemocytometer. A sample of unlysed cells from the same culture was used as a control, from which the lysis percentage was calculated.

# Lysis Rate of *Escherichia coli* with Multiple Passes and Varying Pressure

Lysis of *Escherichia coli* ( $1 \times 10^9$  cells/mL) at Increasing Pressures Over Multiple Passes



**Culture:** A single colony of *Escherichia coli* ATCC 8739 was used to inoculate 10 mL sterile LB (Lysogeny Broth: Tryptone 10 g/L, Sodium Chloride 10 g/L, Yeast Extract 5 g/L) and was grown overnight at 37°C with shaking at approximately 180 rpm. This starter culture was used in a ratio of 1% to inoculate 200 mL sterile LB in 1 L flasks. These were then grown under the same conditions until an Optical Density of approximately 1.2 was reached ( $10^9$  cells/mL).

**Disruption:** The resulting culture was passed in 30 mL aliquots through a Constant Systems Ltd. 'Z+' 1.1 kW Continuous Flow Cell Disruptor at pressures indicated. The machine was rinsed with 30 mL de-ionised water between each use.

**Counting:** After being passed through the machines, 10  $\mu$ L of lysate was mixed 1:1 with the viability stain Trypan Blue. Live cells were counted using a hemocytometer. A sample of unlysed cells from the same culture was used as a control, from which the lysis percentage was calculated.